

medi

Juxta-Fit™ Ankle-Foot Wrap

THIS PACKAGE CONTAINS

(1) Juxta-Fit™ Ankle-Foot Wrap

If there is a problem with this product, please contact CircAid® at (800)-CIRCAID (247-2243) or on the web at www.circaid.com.

Juxta-Fit™ Ankle-Foot Wrap

CARE AND MAINTENANCE



HAND WASH - WARM*



DRIP DRY OR FLAT DRY*



DO NOT BLEACH



DO NOT IRON



DO NOT DRY CLEAN



DO NOT WRING

*To extend the life of the Juxta-Fit™ Ankle-Foot Wrap hand wash and drip or flat dry. The garment can also be machine washed using a gentle cycle and dried on a low heat setting. If doing so, it is recommended that it be washed and dried by itself or in a pillowcase to prevent lint from getting caught in the hook tape. Should lint collect on the hook tape, it can be removed with tweezers.

JF-AFW-DFU-112612-01

JUXTA-FIT™ ANKLE-FOOT WRAP

COMPONENTS OF JUXTA-FIT™ ANKLE-FOOT WRAP:

Polyethylene Plastic, Polyurethane, Lycra®, Nylon, Silver

COMPONENTS OF LABEL:

Polyvinyl Chloride Plastisol or Satin Acetate

THIS PRODUCT CONTAINS SILVER.

THIS PRODUCT DOES NOT CONTAIN LATEX.

INDICATIONS FOR USE: Lymphedema, other forms of Edema: Venous, Post Traumatic, Post Surgical, "Lipedema", Varicose veins (all types), Chronic Venous Insufficiency ("CVI"), Venous Stasis Disease, Venous Valvular Insufficiency ("VVI"), Angio dysplasi, Venous Insufficiency, Post Thrombotic Syndrome, Venous Ulcer (Stasis Ulcer), Post Sclerotherapy, Thrombosis Prevention, Lipodermatosclerosis, Venous Excema.

CAUTIONS: Decreased or absent sensation in the foot, Sensitivity to product components.


CONTRAINDICATIONS: Moderate and/or Severe Peripheral Arterial Disease, Septic Phlebitis, Phlegmasia Cerula Dolens, Decompensated Congestive Heart Failure, allergy to compression materials, infection in the ankle and/or foot.

Custom garments require a prescription. This garment is not designed to be cut or modified, except as outlined in the directions for use; doing other than instructed will void the three month limited warranty. CircAid Medical Products, Inc. assumes no liability for injury or damages resulting from contraindicated use of its products. CircAid is a registered trademark and Juxta-Fit is a trademark of CircAid Medical Products, Inc. Lycra is a registered trademark of INVISTA. ©2012 CircAid Medical Products, Inc. All Rights Reserved.

CIRCAID®
by medi



CircAid Medical Products, Inc.
9323 Chesapeake Drive, Suite B2
San Diego, CA 92123 U.S.A.
01-858-576-3550



MDSS
Schiffgraben 41
30175 Hannover, Germany

DIRECTIONS FOR USE

DONNING THE JUXTA-FIT™ ANKLE-FOOT WRAP:

STEP 1: Don sock or stockinette.

STEP 2: First, unroll the band closest to the toes and secure it at its loosest position to the opposing side of the foot wrap, creating a pocket for the foot to slide through. Slide the foot into the top opening like a slipper, and firmly press the heel into the heel cup.

STEP 3: Ensure that the foot is properly aligned within the garment. The foot piece may need to be pulled toward the front of the foot to insure a proper fit with the heel positioned firmly and snugly into the heel cup.

STEP 4: Tighten the foot-bands on the top of the foot pulling to a firm and comfortable compression.

STEP 5: Detach the long ankle strap and pull it across the front of the foot, wrap with firm and comfortable compression around the ankle, insure there are no gaps, and engage the hook tab to the loop strip. Adjust the ankle strap if any gaps or pressure points occur. **Note: The long ankle strap may pull from a different side depending on which foot, left or right, the AFW is applied to.**

STEP 6: Readjust the straps as necessary throughout the day.

Note: Slightly loosen the bands for night-time wear. If the standard length of the foot piece is too long for a proper fit please see the next page for foot length adjustment instructions.

REMOVAL: Unwrap the long ankle-band, pass it underneath the arch of the foot and reattach it to the body of the garment. Loosen the foot-bands enough so that the foot can easily slide out of the garment. When complete, the garment should resemble its starting position for easy reapplication.

AFW READY FOR REAPPLICATION



STEP 2



STEP 3



STEP 4



STEP 5



FINISHED

ANKLE-FOOT WRAP (AFW) FOOT-LENGTH ADJUSTMENT:

The following steps will permanently reduce the AFW's foot-length by 3 cm (1.25").


STEP 1: Locate the two notches on both sides of the Juxta-Fit™ Ankle-Foot Wrap approximately 3 cm (1.25") from the front of the garment.

STEP 2: Draw a straight line from one notch to the other.

STEP 3: Carefully cut along the drawn line, assuring that the reinforced seam section of the hook material is not cut.

PLEASE NOTE: NO OTHER MODIFICATIONS CAN BE MADE TO THE JUXTA-FIT™ ANKLE-FOOT WRAP.

BE CAREFUL THAT THE REINFORCED SEAM WITHIN THE HOOK MATERIAL IS NOT CUT WHILE SHORTENING THE GARMENT.



REMEMBER: The Juxta-Fit™ Ankle-Foot Wrap should never hurt. It should feel firm but comfortable. If pressure increases or decreases noticeably during wear, loosen or tighten any bands accordingly.

IF EXPERIENCING ANY PAIN, IMMEDIATELY REMOVE THE JUXTA-FIT™ ANKLE-FOOT WRAP AND CONSULT YOUR PHYSICIAN.

CIRCAID
by medi

If further assistance is required, contact CircAid® Customer Service by phone at (800) CIRCAID (247-2243) or by email at info@circaid.com