JUXTA FITTM

ARMSLEEVE

TM

medi

BUILT-IN PRESSURE SYSTEMTM (BPSTM)

Pull bands to firm and comfortable compression, or if directed by a Healthcare Professional, follow steps below to set compression to the prescribed level.

USING THE BUILT-IN PRESSURE SYSTEMTM (BPSTM):

STEP 1: Locate the BPSTM card in your packaging.

STEP 2: Starting with the bottom band, line up the pressure guide line on the garment with the line on card marked with ▼.


STEP 3: Note where the second pressure guide line on the garment is in relation to the correct color coded pressure range on the card. (Fig. A, B, and C)

STEP 4: If the pressure guide line on the garment either falls short (Fig. B) or goes beyond the correct pressure range (Fig. C), readjust the band as necessary so the second guideline is aligned with the prescribed pressure range (Fig. A). Completely loosen band before reapplying if too much tension has been applied.

STEP 5: Repeat BPSTM Steps 2-4 for each band going up the arm. Be sure to adjust to the SAME tension on the BPSTM card scale as the bottom band. Bands may need adjusting throughout the day to maintain the prescribed pressure range.

This garment is designed to provide measurable gradient compression based on average compression levels at the wrist.

CircAid Medical Products, Inc. assumes no liability for injury or damages resulting from contraindicated use of its products. CircAid is a registered trademark and BPS, Built-In Pressure System, and Juxta-Fit are trademarks of CircAid Medical Products, Inc. ©2012 CircAid Medical Products, Inc. All Rights Reserved.

▲ Fig. A Correctly Adjusted for 20-30mmHg

▲ Fig. B Not enough tension for 20-30mmHg

▲ Fig. C Too much tension for 20-30mmHg, Set to 40mmHg+

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