# JUXTA FIT<sup>TM</sup>

## ESSENTIALS

### UPPER LEG WITH KNEE

TM

medi

## BUILT-IN PRESSURE SYSTEM<sup>TM</sup> (BPS<sup>TM</sup>)

Pull bands to firm and comfortable compression, or if directed by a Healthcare Professional, follow steps below to set compression to the prescribed level.

### USING THE BUILT-IN PRESSURE SYSTEM<sup>TM</sup> (BPS<sup>TM</sup>):

**STEP 1:** Locate the Built-In Pressure System<sup>TM</sup> guide card. Highlight the edge of the card that correctly references the patient's ankle circumference measurement.

**STEP 2:** Starting with the bottom band, line up the pressure guide line on the garment with the line on card marked with ▼.

**STEP 3:** Note where the second pressure guide line on the garment is in relation to the correct color coded pressure range on the card. (Fig. A, B, and C)

**STEP 4:** If the pressure guide line on the garment either falls short (Fig. B) or goes beyond the correct pressure range (Fig. C), readjust the band as necessary so the second guideline is aligned with the prescribed pressure range (Fig. A). Completely loosen band before reapplying if too much tension has been applied.


**STEP 5:** Repeat BPS<sup>TM</sup> Steps 2-4 for each band going up the leg. Be sure to adjust to the SAME tension on the BPS<sup>TM</sup> card scale as the bottom band. Bands may need adjusting throughout the day to maintain the prescribed pressure range.

\*This garment is designed to provide gradient compression based on average ankle compression levels.

CircAid Medical Products, Inc. assumes no liability for injury or damages resulting from contraindicated use of its products. CircAid is a registered trademark and BPS, Built-In Pressure System, and Juxta-Fit are trademarks of CircAid Medical Products, Inc. ©2012 CircAid Medical Products, Inc. All Rights Reserved.

▲ Fig. A Correctly Adjusted for 30-50mmHg

▲ Fig. B Not enough tension for 30-50mmHg

▲ Fig. C Too much tension for 30-50mmHg, Set to 50mmHg+

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