

medi

## Juxta-Fit™ Interlocking Ankle-Foot Wrap

THIS PACKAGE CONTAINS

One (1) Juxta-Fit™ Interlocking Ankle-Foot Wrap

## Juxta-Fit™ Interlocking Ankle-Foot Wrap

CARE AND MAINTENANCE



HAND WASH – WARM\*



DRIP DRY OR FLAT DRY\*



DO NOT BLEACH



DO NOT IRON



DO NOT DRY CLEAN



DO NOT WRING

\*To extend the life of the Juxta-Fit™ Interlocking Ankle-Foot Wrap hand wash and drip or flat dry. The garment can also be machine washed using a gentle cycle and dried on a low heat setting. If doing so, it is recommended that it be washed and dried by itself or in a pillowcase to prevent lint from getting caught in the hook tape. Should lint collect on the hook tape, it can be removed with tweezers.

JF-IAFW-DFU-112612-01

# JUXTA-FIT™

## INTERLOCKING ANKLE-FOOT WRAP

**COMPONENTS OF JUXTA-FIT™ INTERLOCKING ANKLE-FOOT WRAP:**

Polyethylene Plastic, Polyurethane, Lycra®, Nylon, Silver

**COMPONENTS OF LABEL:** Polyvinyl Chloride Plastisol or Satin Acetate

**THIS PRODUCT CONTAINS SILVER.**

**THIS PRODUCT DOES NOT CONTAIN LATEX.**

**INDICATIONS FOR USE:** Lymphedema, other forms of Edema: Venous, Post Traumatic, Post Surgical, "Lipedema", Varicose veins (all types), Chronic Venous Insufficiency ("CVI"), Venous Stasis Disease, Venous Valvular Insufficiency ("VVI"), Angio dysplasi, Venous Insufficiency, Post Thrombotic Syndrome, Venous Ulcer (Stasis Ulcer), Post Sclerotherapy, Thrombosis Prevention, Lipodermatosclerosis, Venous Excema.

**CAUTIONS:** Decreased or absent sensation in the foot, Sensitivity to product components.

**CONTRAINDICATIONS:** Moderate and/or Severe Peripheral Arterial Disease, Septic Phlebitis, Phlegmasia Cerula Dolens, Decompensated Congestive Heart Failure, allergy to compression materials, infection in the ankle and/or foot.

This garment is not designed to be cut or modified in any way, doing so will void the three month limited warranty. CircAid Medical Products, Inc. assumes no liability for injury or damages resulting from contraindicated use of its products. CircAid is a registered trademark and Juxta-Fit is a trademark of CircAid Medical Products, Inc. Lycra is a registered trademark of INVISTA. Velcro is a registered trademark of Velcro USA, Inc. ©2012 CircAid Medical Products, Inc. All Rights Reserved.

CIRCAID®  
by medi



CircAid Medical Products, Inc.  
9323 Chesapeake Drive, Suite B2  
San Diego, CA 92123 U.S.A.  
01-858-576-3550



MDSS  
Schiffgraben 41  
30175 Hannover, Germany

# DIRECTIONS FOR USE

## DONNING THE JUXTA-FIT™ INTERLOCKING ANKLE-FOOT WRAP (AFW):

**STEP 1:** Don sock or stockinette.

**STEP 2:** First, unroll the band closest to the toes and secure it at its loosest position to the opposing side of the foot wrap. Ensure that the interlocking heel-bands are secured at their loosest position. Slide the foot into the top opening like a slipper, and firmly press the heel into the heel opening.

**STEP 3:** Ensure that the foot is properly aligned within the garment; the front edge of the foot wrap should align with the base of the toes. Adjust the interlocking heel-bands to fit snugly around the back of the heel.

**STEP 4:** Tighten the foot-bands on the top of the foot, pulling to a firm and comfortable compression.

**STEP 5:** Detach the long ankle-band, pull it across the front of the foot, wrap with firm and comfortable compression around the ankle, ensure there are no gaps, and engage the hook tab to the loop strip. Adjust the ankle-band if any gaps or pressure points occur.



**DO NOT TUCK BANDS. THE VELCRO® HOOKS SHOULD NEVER TOUCH THE SKIN!**




**Note:** The long ankle-band may pull from a different side depending on which foot, left or right, the AFW is applied to. Make sure to secure the long ankle-band last, over the legging, if using the AFW with a Juxta-Fit™ legging.



**STEP 7:** Readjust the bands as necessary throughout the day to maintain a firm and comfortable compression level. **Note:** Slightly loosen the bands for night-time wear.

**REMOVAL:** Unwrap the long ankle-band, pass it underneath the arch of the foot and reattach it to the body of the garment. Loosen the bands on top of the foot enough so that the foot can easily slide out of the garment. When complete, the garment should resemble its starting position for easy reapplication.



**REMEMBER:** The Juxta-Fit™ Interlocking Ankle-Foot Wrap (AFW) should never hurt. It should feel firm but comfortable. If pressure increases or decreases noticeably during wear, loosen or tighten any bands accordingly.



**IF EXPERIENCING ANY PAIN, IMMEDIATELY REMOVE THE JUXTA-FIT™ INTERLOCKING AFW AND CONSULT YOUR PHYSICIAN.**