

# COMFORT™

## EZ SINGLE-BAND ANKLE-FOOT WRAP

**COMPONENTS OF COMFORT EZ SINGLE-BAND ANKLE-FOOT WRAP™:** Nylon and Spandex

**THIS PRODUCT DOES NOT CONTAIN LATEX.**

**INDICATIONS FOR USE:** Lymphedema, other forms of Edema: Venous, Post Traumatic, Post Surgical, "Lipedema", Varicose veins (all types), Chronic Venous Insufficiency ("CVI"), Venous Stasis Disease, Venous Valvular Insufficiency ("VVI"), Angio dysplasi, Venous Insufficiency, Post Thrombotic Syndrome, Venous Ulcer (Stasis Ulcer), Post Sclerotherapy, Thrombosis Prevention, Venous Excema, Lipodermatosclerosis.

**CAUTIONS:** Decreased or absent sensation in the foot, Sensitivity to product components.

**CONTRAINDICATIONS:** Moderate and/or Severe Peripheral Arterial Disease, Septic Phlebitis, Decompensated Congestive Heart Failure, Phlegmasia Cerula Dolens, Allergy to compression materials, Infection in the ankle and/or foot.

CircAid Medical Products, Inc. assumes no liability for injury or damages resulting from contraindicated use of its products. CircAid is a registered trademark and Comfort EZ Single-Band Ankle-Foot Wrap is a trademark of CircAid Medical Products, Inc.

©2012 CircAid Medical Products, Inc. All Rights Reserved.



CircAid Medical Products, Inc.  
9323 Chesapeake Drive, Suite B2  
San Diego, CA 92123 U.S.A.  
01-858-576-3550

**CIRCAID**  
by medi

EC REP

MDSS  
Schiffgraben 41  
30175 Hannover, Germany

C-EZAFW-DFU-111512-01

## DIRECTIONS FOR USE:

**BEFORE DONNING THE COMFORT EZ SINGLE-BAND ANKLE-FOOT WRAP™:** Put on a sock. The Comfort EZ Single-Band Ankle-Foot Wrap™ (AFW) should not be applied directly against the skin

**DONNING THE COMFORT EZ SINGLE-BAND ANKLE-FOOT WRAP™:**

**STEP 1.** Place the hook end of the AFW on top of the foot, facing away from the skin. Holding the hook tab in place, wrap the AFW around the foot once, and secure the AFW to the hook end.



STEP 1

**STEP 2.** Continue wrapping the AFW around the foot with firm and comfortable compression. Cover the foot from the base of the toes to the heel making sure there are no gaps between the layers of elastic. Apply double layers as desired to target problem areas.



STEP 2

**STEP 3.** If applying for the first time, the band can be shortened as needed. Secure the loose end of the AFW using one of the provided hook tabs.

*Remember: The Comfort EZ Single-Band Ankle-Foot Wrap™ should never hurt. It should feel firm but comfortable. If pressure increases or decreases noticeably during wear, readjust the AFW as needed.*

**REMOVING THE COMFORT EZ SINGLE-BAND ANKLE-FOOT WRAP™:**

Disengage the hook tab and remove it from the AFW. Unwrap the foot and roll the AFW into a cylinder to prepare for the next application. Make sure the sewn-in hook tab is at the tail-end of the roll.



STEP 3

## CARE INSTRUCTIONS:



**MACHINE WASH WARM\***



**TUMBLE DRY LOW\***



**DO NOT BLEACH**



**DO NOT IRON**



**DO NOT DRY CLEAN**



**DO NOT WRING**



FINISHED

The Comfort EZ Single-Band Ankle-Foot Wrap™ can be machine washed with mild detergent using a gentle cycle and machine dried on permanent press or a low heat setting. It is recommended that it be washed and dried by itself or in a pillowcase to prevent lint from getting caught in the hook end. To extend the life of the Comfort EZ Single-Band Ankle-Foot Wrap™ hand wash and drip or flat dry.