

CARE AND MAINTENANCE



HAND WASH



DRIP OR LINE DRY



DO NOT MACHINE WASH



DO NOT BLEACH



DO NOT IRON



DO NOT DRY CLEAN



DO NOT WRING

The Comfort Non-Slip Liner™ should be hand washed and air dried. Due to the nature of the material, debris will adhere to it; use tweezers to remove debris that does not wash out. Note that debris will decrease the effectiveness of the material over time.



This product should not be used by children six years of age and younger.



COMFORT NON-SLIP LINER™

The Comfort Non-Slip Liner™ is worn directly against the skin, underneath any compression legging. Its tacky surface aids in preventing the garment from sliding.

COMPONENTS OF COMFORT NON-SLIP LINER™:
Polyester fiber with binders.


THIS PRODUCT DOES NOT CONTAIN LATEX.

CircAid Medical Products, Inc. assumes no liability for injury or damages resulting from contraindicated use of its products. CircAid is a registered trademark. Comfort Non-Slip Liner is a trademark of CircAid Medical Products, Inc. ©2012 CircAid Medical Products, Inc. All Rights Reserved.

CIRCAID[®]
by medi



CircAid Medical Products, Inc.
9323 Chesapeake Drive, Suite B2
San Diego, CA 92123 U.S.A.
01-858-576-3550



MDSS
Schiffgraben 41
30175 Hannover, Germany

DIRECTIONS FOR USE

TRIMMING THE COMFORT NON-SLIP LINER™

(For use with your Juxta-Fit™ Legging)

STEP 1 Unfold the Comfort Non-Slip Liner™ and place it flat on a surface. *Be sure to choose a clean surface as the material attracts debris.*





Unroll the tab ends of the legging. Place the legging, or portion of the legging (typically the thigh portion), flat on top of the Comfort Non-Slip Liner™ with the black side of the material facing up.

Note: Do not let the hook tabs come in contact with the Comfort Non-Slip Liner™. The hook tabs may damage the Comfort Non-Slip Liner™ and the fibers of the Comfort Non-Slip Liner™ will decrease the effectiveness of the hook tabs.

STEP 2 Using the legging as a template, cut a straight line through the Comfort Non-Slip Liner™. Use the outside edge of the tab ends as a guide.

BE SURE NOT TO CUT THE LEGGING, AS THIS WILL AUTOMATICALLY VOID THE LEGGING'S WARRANTY.

STEP 3 If needed, trim the Comfort Non-Slip Liner™ along the top and bottom edge of the legging to minimize excess material.



TRIMMING THE COMFORT NON-SLIP LINER™


(Using Circumference Measurements)

STEP 1 Measure the height of the desired coverage area on the limb.

Note: The maximum height of coverage area is 46cm.

STEP 2 Measure the largest and bottom circumference of the limb in the coverage area. Add 25cm to each of these measurements to ensure your Comfort Non-Slip Liner™ overlaps.

STEP 3 Using these three measurements, mark and cut the Comfort Non-Slip Liner™ per the drawing below.



COMFORT NON-SLIP LINER™ WITH JUXTA-FIT™ LEGGING

DONNING THE COMFORT NON-SLIP LINER™

STEP 1 Wrap the Comfort Non-Slip Liner™ around the part of the limb where the legging will be donned (e.g. thigh). The Comfort Non-Slip Liner™ should be worn directly against skin; it will not be effective if worn over sock liners.

STEP 2 Overlap the excess material at the front of the limb to anchor the Comfort Non-Slip Liner™ in place. The Comfort Non-Slip Liner™ will adhere to itself; no tape or other adhesives are necessary.

STEP 3 Don the legging over the Comfort Non-Slip Liner™. Be sure that none of the hook tabs are adhered to the Comfort Non-Slip Liner™. The hook tabs should be entirely secured to the garment. Getting the hook tabs tangled in the Comfort Non-Slip Liner™ may damage the Comfort Non-Slip Liner™ as well as decrease the effectiveness of the hook tabs.



If further assistance is required, please contact CircAid® Customer Service by phone at (800) CIRCAID (247-2243) or by email at info@circaid.com.